



**MODEL NO. - ACM- 82302**

# WILLY MILL

# WILLY MILL

## INTERMEDIATE MODEL 'ARTHUR THOMAS TYPE'

Electrically operated to work on 220 Volts A.C Recommended for grinding of dry plant samples. Also for a great variety of commercial dry materials. Unit having two stationary steel blades and Rotary with steel edges provided in the cast Iron chamber. A removable safety Glass sheet plate from the face of chamber of permits observation during operation. Willey Mill is fitted on substantial Cast iron base. Complete unit supplied with aluminum / Cast chamber , polished hopper with cover, wooden plugger, Three delivery tubes mounted with sieves cloth mesh no; 20, 40, and 60mesh , One Allen Key for adjusting blades and 3 Glass Jars, brush etc.



## **WORKING CHAMBER**

- 40 mm x 20 mm
- 65 mm x 25 mm
- 100 mm x 50 mm

## **WILLY MILL (HEAVY DUTY)**

WITH 1 H.P. Electric Motor to work on 220 Volt A.C

Extra Heavy Duty mill mounted in enclosed iron frame. Six Bolted with grinding chamber. A set of 3 Brass Sieve No; 2mm , 4mm, 6 mm .Complete unit supplied with Heavy floor covered with Iron sheet.



An ISO 9001 : 2208 | ISO 14001 : 2008 | ISO 13485  
WHO: GMP Products | GLP Compliant Products

[www.acmasindia.com](http://www.acmasindia.com) | [www.acmasglobal.com](http://www.acmasglobal.com) | [www.test-chambers.com](http://www.test-chambers.com)

## **ACMAS TECHNOLOGIES PVT. LTD.**

### **CORPORATE OFFICE/ SHOWROOM (INDIA)**

Plot No. 352-353, Sector – 57  
Phase- IV Kundli, Sonapat, Haryana - 131028  
Land Line No.: +91 - 11- 47619688,  
Fax: +91-11-47619788  
E-mail.: [info@acmasindia.com](mailto:info@acmasindia.com)

### **SALES OFFICE (HONG KONG)**

Unit D 28 11/F Wing Tat Comm,Bidg 97, Bonham Strand East,  
Sheung Wan, Hong Kong (PRC)  
Tel.: 0086-13929598046 | 0086- 18922303099  
E-mail.: [hk@acmasindia.com](mailto:hk@acmasindia.com)

### **SALES OFFICE (RUSSIA)**

Inmed Trade Street Ozerkovsky Embankment,  
Unit No 50, Straine- 1, Off- 502, Moscow, Russia  
E-mail.: [info@acmasindia.com](mailto:info@acmasindia.com)  
Tel.: 0049- 79592345 | Email: [russia@acmasindia.com](mailto:russia@acmasindia.com)